



- Shear-free micro wells lead uniform & perfect spheroid cell culture
- Easy to handle & cutting down on labor intensiveness
- Individual sterilized package





### Shear-free • Identical Size Control • High-throughput




#### SPECIFICATION


Area of wells	29 x 29mm
Well size	500um(I.D.) x 500um(Depth)
Color	Clear
No. of wells	361 micro-wells
Model	ISF-361
Material	PDMS(Polydimethylsiloxane)



#### Conventional 3D cell culture plate



#### INCYTO SpheroFilm



#### ADVANTAGES

1. Extensible culture period of 10-30 days due to easy access of medium change.
2. Extremely shortened culture period of 2-3 days for some application by controlling cell seeding density.
3. Co-culture ready system without direct cell-cell contact.
4. Generating size-controlled cell spheroids.
5. Easy to change the growth medium in the device due to shear stress-free structure.
6. Applicable to various cell types.
7. Unlimited throughput (361 cell spheroids can be obtained from one device).

